



1

Compressive Strength

◆ Basic Concept

Amalgam is strong in compression

Weak in tension & shear → (Brittle material)



One-Hour Compressive Strength (According to ANSI/ADA)

Alloy Type	1-hour Compressive Strength
Lathe-cut Low Copper	45 MPa <small>مش مهم الأرقام</small>
Spherical Low Copper	88 MPa
Admixed High Copper	118 – 140 MPa
High Copper Unicompositional	250 MPa
ANSI/ADA Requirement	AT LEAST 80 MPa

Compressive strength

- **1 hour** If strength < 80 MPa →
 - ⚠ Risk of:
 - Cracks
 - Fracture
 - Failure due to premature biting
- **7 days** Highest compressive strength = Unicompositional high copper alloys

Tensile Strength

◆ Key Facts

- Tensile strength = **only a fraction** of compressive strength
- Unicompositional high copper alloy → **Highest tensile strength**
- Indicates strength of amalgam phases

Strong Phase Order: $\gamma > \gamma_1 > \gamma_2$

- γ (Unreacted Ag_3Sn) = Strongest
- γ_1 (Ag_2Hg_3) = Intermediate
- γ_2 (Sn_8Hg) = Weakest & most corrodible



⚠ Effect of Mercury Amount on Amalgam Phases

💡 If Excess Mercury (Hg) is Used:

- Excess Hg → ↑ Formation of γ_1 and γ_2 phases
- This occurs at the expense of the unreacted γ phase

💡 What Happens Inside the Mix?

- More γ_1 (Ag_2Hg_3) forms
- More γ_2 (Sn_8Hg) forms
- ↓ Amount of unreacted γ (Ag_3Sn) phase

💡 Relation to Strength of Amalgam

☒ With Excess Hg:

- $\rightarrow \gamma_1 + \gamma_2 \uparrow$
- \rightarrow Unreacted $\gamma \downarrow$
- ✗ Result: Weaker amalgam mass

🏆 With Proper (Optimal) Hg Amount:

- \rightarrow Less γ_1 and γ_2 formation
- \rightarrow More unreacted γ phase remains
- ◆ Result: Stronger amalgam mass

Elastic Modulus (Stiffness)

- Definition : 🧠 Indicates stiffness of amalgam.

📌 Facts

- High copper alloys → Higher modulus of elasticity ($\approx 11-20$)
- Therefore:
- 👉 High copper amalgam = More rigid & stiffer than low copper

Component	Elastic Modulus (GPa)
Dentin	14-18.6
Enamel	80

Creep (Viscoelastic Behavior)

- Definition

Permanent deformation under static load below elastic limit.

Creep Comparison

ANSI/ADA Maximum Allowable Creep = 3%

Alloy Type	Creep Value
Lathe-cut Low Copper	Highest creep ✗
Unicompositional High Copper	Least creep ☑

🎯 Factors Affecting Creep

- ↓ Number of η crystals
- ↑ Mercury content
- Grain size of γ_1 phase
- Volume of γ phase



Dimensional Changes

Pattern

- Initial contraction بسب Hg dissolution → then expansion بسب Crystal formation in different amalgam phases

Timeline

- Stops after: 6–8 hours
- Final value reached: 24 hours
- Final change usually: Negative (slight contraction)

Best Alloy

◆ High Copper (Unicompositional) → Least dimensional change

◆ ADA Permissible Dimensional Change: $\pm 20 \mu\text{m}/\text{cm}$

Excessive Delayed Expansion (VERY IMPORTANT EXAM POINT)

◆ Occurs When: ➔ Zinc-containing alloy + Moisture contamination

◆ Chemical Reaction: $\text{H}_2\text{O} + \text{Zn} \rightarrow \text{ZnO} + \text{H}_2$ (gas)

◆ Consequences:

- Blistering of amalgam surface
- Excessive delayed expansion
- Pain (pressure $\approx 400 \mu\text{m}/\text{cm}$)
- Corrosion later

⌚ Reaction occurs: 24–72 hours after insertion

Clinical Tip

If moisture contamination is expected (e.g., Children):  Use Zinc-free alloy

⚠ Even in zinc-free:

Moisture contamination:

- Prevents cohesion of amalgam layers
- Leads to formation of a **weak lamellated mass**



Corrosion of Amalgam

- ◆ Definition : Destruction of metal by chemical/electrochemical reaction

Effects of Corrosion

- 1 Increased porosity
- 2 Reduced marginal integrity
- 3 Loss of strength
- 4 Metallic products lead to a metallic taste in the mouth

Corrosion of Phases

- Least corrodible phase → γ_1
- Most corrodible phase → γ_2 (Sn_8Hg)

Feature	▼ Corrosion in Low Copper	▲ High Copper Amalgam
Main Site of Corrosion	Occurs throughout γ_2 phase (network structure)	No γ_2 phase present
Phases Affected	γ_2 (Sn_8Hg) network phase	High copper phases: η & eutectic Ag-Cu
Reaction Products	Tin oxychloride + Mercury compounds	$\text{CuCl}_2 \cdot 3\text{Cu}(\text{OH})_2$
Structural Nature	γ_2 is a network structure → more corrosion	η phase is NOT network like γ_2
Effect on Strength	Causes porosity & ↓ strength of set amalgam	Higher resistance to corrosion → better strength
Overall Corrosion Resistance	Low ✗	High ☑

Manipulation of Dental Amalgam (Step by Step)


1 Selection of Alloy

- 🔍 Factors Affecting Selection • Particle size :


Coarse cut	Fine cut	Microcut
(750 μm) Difficult to dispense and difficult to carve, it leaves rough surface.	(15–35 μm) easy to dispense, better on carving, leaves smooth surface.	(less than 3 μm) very large surface area per volume needs too much Hg to wet all particles → ↑ Hg ratio leads to ↑ setting time, ↓ strength, ↑ creep, ↑ setting expansion, ↑ tarnish and corrosion.



Particle Shape

 Lathe Cut Particles	 Spherical Particles
Irregular shape with large surface area per volume → needs more Hg for mixing	Lower surface area → needs less Hg for mixing
Offers high resistance to condensation forces	Less resistance to condensation forces
High friction between particles → better marginal integrity	Less friction between particles → less marginal integrity
Rough surface requires more	after carving → finishing
	

Spheroidal shaped : cucumber gives good marginal Integrity



Particle Composition

According to Copper Content:

- Low Copper Alloy
- High Copper Alloy (BEST)

Advantages of High Copper:

- Eliminates γ_2 phase
- High early strength
- Low creep
- Better corrosion resistance
- Less marginal fracture

Type	Zinc %
Zinc-containing	> 0.01% Zn
Zinc-free	< 0.01% Zn

According to Zinc Content

 Zinc alloy must NOT be contaminated with moisture→ To avoid delayed expansion



2 Proportioning of Mercury to Alloy

Must be Accurate Because:

X Excess Mercury

- ↑ γ_1 & γ_2 formation
- ↑ Setting time
- ↓ Strength
- ↑ Creep
- ↑ Tarnish & Corrosion

X Less Mercury

- Poor wetting of particles
- Friable mix
- More voids
- ↓ Strength
- ↑ Corrosion liability





Ideal Mercury Ratio



- Lathe-cut particles → 54% Hg
 - Spherical particles → 43% Hg
- 👉 preweighed tablets of 400 mg is used to facilitate proportioning (used with hand trituration)
- 👉 preweighed proportioned capsules that contain proper amount of alloy and mercury are used with mechanical Amalgamators .

3 Trituration (Mixing)

Definition

Process where mercury reacts with alloy powder to form a **homogenous, shiny, plastic mass**.

Types of Trituration

Hand Trituration



Mechanical Trituration

- **Amalgamizer** : which provides Automatic dispensing and mechanical mixing .
- **Amalgamator** : which mixes proportioned capsules , in which the capsule is attached to capsule holder , which is attached to a motor that allows eccentric movements giving the force essential for mixing, 7 10 s.

Objectives of Trituration

- 1 Remove oxide film from alloy particles
- 2 Coat particles with mercury
- 3 Supply energy for amalgamation reaction

Outcomes of Amalgamation Reaction

- ↓ Unreacted γ phase size
- Formation of γ₁ phase (non-corrodible)
- Formation of γ₂ phase (in conventional alloy)
- Progressive Hg consumption → ↓ plasticity



Mulling Operations

- It aims to produce more homogenous mix with good texture.
- **Hand trituration:** by gentle rubbing of the mix between fingers.
- **Mechanical trituration:** by further amalgamation 2–3 seconds at low speed.

⚠ Over vs Under Trituration

🔥 Over-trituration

- ↓ Working time
- ↑ Creep
- Slight ↑ contraction

❄️ Under-trituration

- More voids
- ↓ Strength
- ↑ Corrosion liability



Condensation (VERY CLINICAL)

🧠 **Definition :** Packing amalgam into cavity **layer by layer** using condensers.

- During this step the material is highly adapted to the cavity walls

⌚ Time required ≈ 3 minutes

- the amount of mercury in the final restoration is determined in this stage



amalgam carrier carrying amalgam into a prepared cavity layer by layer



Condensation of amalgam; a: small condenser, followed by exceedingly larger condenser (b) to overpack restoration

◆ Principal of condensation:

Condense the amalgam mass and adapt it to the **preparation walls, floor, and matrix** (when used).

◆ Purpose:

- Produces a restoration free of voids → ↑ density & ↓ marginal leakage
- Minimizes mercury content → ↓ corrosion & ↑ restoration strength
- Ensures cohesion with subsequent increments & prevents layering
- Facilitates blotting of excess mercury for better final restoration properties



Methods of Condensation

Hand Condensation

- Adjusted according to cavity size and particle size
- Lathe cut alloys: small diameter (1–2 mm) + large vertical force
- Spherical alloys: large diameter + less lateral force

Mechanical Condensation

- Based on impacts or vibration
- More effective for lathe cut alloys than spherical alloys
-  Not recommended for spherical alloys → can spread droplets & ↑ mercury vapor

Effects of Delayed Condensation

- Amalgam partially sets → loses plasticity → poor adaptation to cavity walls
- Excess mercury cannot be removed →
- ↓ compressive strength
- ↑ corrosion susceptibility
- ↑ setting expansion
- ↑ creep

Moisture Contamination

- Zinc-containing alloys: release H₂ → delayed expansion
- Zinc-free alloys: weakened by lamellation

5 Pre-Carve Burnishing

- ◆ Tool: Large ball burnisher
- ◆ Purpose:
- ↑ adaptation of condensed amalgam to cavity margins
- Facilitates carving
-  Objection: may carry excess Hg into critical CSA



6 Carving

- ◆ Definition: Procedure to reproduce anatomical & functional relations of restored tooth


Timing & Technique

- After initial setting
- Direction: Any direction except from restoration to tooth → prevents submargins
- Carver: Must be sharp → dull carver induces stress, cracks & Hg distribution issues





Steps



1. **Initial carving:** Use carving explorer to reproduce occlusal embrasures
2. **Final carving:** Sharp carver parallel to cusp inclined planes using **discoid, cleoid, diamond-shaped, or Hollenback carvers**

7 Post-Carve Burnishing

- **Aim:** Smooth restoration, enhance adaptation
- **Technique:** Light rubbing with anatomical burnisher from **center → cavity margins**, follow anatomical landmarks
- **Appearance:** Satin, not shiny
- **Timing:** After carving, before finishing & polishing



8 Finishing & Polishing

Timing

- Low copper alloys → after 24 hours
- High copper alloys → 8–10 min (high early strength)

Objectives

- Avoid rough surfaces → prevent micropits → ↓ galvanic action → ↓ tarnish & corrosion
- Remove marginal flashes/overhangs, correct contacts & contours
- Produce **lustrous, homogeneous surface** → ↑ corrosion resistance, ↓ plaque retention, ↑ biocompatibility, ↓ gingival irritation, ↑ stress response
- Encourage better home care

Limited Lifespan of Amalgam Restoration

Survival Factors

- Inversely proportional to restoration size
- Failures mainly due to:
 - Creep
 - Corrosion
 - Marginal breakdown

◆ Preferred Choice in Modern Dentistry:
→ **High Copper Amalgam** (Low creep + High corrosion resistance)





Mercury Safety & Hygiene (Clinical Must-Know)

- Ensure good ventilation
 - Use sealed capsules
 - Clean spills immediately (NO vacuum cleaner)
 - Wear gloves, mask, eye protection
 - Proper waste disposal
 - Avoid direct skin contact with amalgam mix
-



Instruments Required for Amalgam Restoration

- Diagnostic set (Mirror, Tweezer, Explorer)
 - Amalgam carrier
 - Condensers (small & large)
 - Burnishers
 - Carvers
 - Mortar & pestle / Amalgamator
 - Infection control kit (Gloves, Mask, Napkin)
- 